

LONGER RS485 Protocol for BT100-1F

1. Data Format: 1start + 8data + 1even parity + 1stop, 1200bps.

This defines the data format: 1 start bit, 8 data bits, one even parity bit, and one stop bit at 1200 bits per second.

2. Command Format: flag+ addr + len + pdu + fcs.

flag: E9H is the start **flag** of a command string. Every command string is preceded with the start of E9H.

- In one command string, there is no other E9H except start **flag** E9H. When transmitting, E8H is replaced by E8H 00H, and E9H is replaced by E8H 01H except start **flag**. When receiving, E8H 00H is replaced by E8H, and E8H 01H is replaced by E9H.

addr: Pump address (i.e. Pump I.D.#.), take up 1 byte.

- The pump address can be set from 1 to 30. 31(1F) is broadcast address.
- In a command string from the control computer, if the **addr** is pump address, the corresponding pump will execute the command and respond. And if the **addr** is broadcast address, all the pumps execute the same command, and pumps don't respond.


len: Length of **pdu**, take up 1 byte.

Fcs: XOR of **addr**, **len**, **pdu**, take up 1 byte.


3. Pdu Format: application layer code format

3.1. Write Running Parameter of Flow Control Mode

Control computer command string:



Pump response:



- Flow unit: nL, 1 L=10³ mL=10⁶ μL=10⁹ nL.

- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31). When the **addr** is pump address, the corresponding pump will execute the command and respond. When the **addr** is broadcast address, all the pumps execute the same command, and pumps don't respond.

LONGER RS485 Protocol for BT100-1F

3.2. Read Running Parameter of Flow Control Mode

Control computer command string:

RF

Pump response:

RF Show flow (4 bytes) State1 (1byte)

- Flow unit: nL, 1 L=10³ mL=10⁶ μL=10⁹ nL.
- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

3.3. Write Dispensing Parameter

Control computer command string:

WD Disp Vol (4 bytes) Copy No. (2 bytes) Disp Flow (4 bytes) Pause time (2 bytes)

Pump response:

WD

- Disp Vol unit: 0.01mL (1) Disp Vol range: 0.01-9990 mL (1-999000)
- Copy No. range: 0-9999, 0 for infinity
- Disp Flow unit: 1nL/min (1) Disp Flow range: 1-1000000000
- Pause time unit: 0.1s (1) Pause time range: 0-59940
- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31). When the **addr** is pump address, the corresponding pump will execute the command and respond. When the **addr** is broadcast address, all the pumps execute the same command, and pumps don't respond.

3.4. Read Dispensing Parameter

Control computer command string:

RD

Pump response:

RD Disp Vol (4 bytes) Copy No. (2 bytes) Disp Flow (4 bytes) Pause time (2 bytes)

- Disp Vol unit: 0.01mL (1) Disp Vol range: 0.01-9990 mL (1-999000)
- Copy No. range: 0-9999, 0 for infinity
- Disp Flow unit: 1nL/min (1) Disp Flow range: 1-1000000000
- Pause time unit: 0.1s (1) Pause time range: 0-59940
- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

LONGER RS485 Protocol for BT100-1F

3.5. Write Running Parameter of Dispensing Control Mode

Control computer command string:



Pump response:



- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31).

3.6. Read Running Parameter of Dispensing Control Mode

Control computer command string:




Pump response:




- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

3.7. Write Pump Head and Tubing.

Control computer command string:



Pump response:




- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31).

3.8. Read Pump Head and Tubing.

Control computer command string:



Pump response:




- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

LONGER RS485 Protocol for BT100-1F

3.9. Write Back Suction Parameter

Control computer command string:




Pump response:




- Back suction time unit: 0.1s (1) Back suction time range: 0-99.9s (0-999)
- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31).

3.10 Read Back Suction Parameter

Control computer command string:




Pump response:




- Back suction time unit: 0.1s (1) Back suction time range: 0-99.9s (0-999)
- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

3.11. Write Pump Address

Control computer command string:



Pump response:



- In a command string from the control computer, the **addr** can be pump address (1-30) and broadcast address (31).
- Pump address can be set one by one with broadcast address.

3.12. Read Pump Address

Control computer command string:



Pump response:



- In a command string from the control computer, if the **addr** is one pump's address (1-30), the corresponding pump will respond.

LONGER RS485 Protocol for BT100-1F

APPENDIX

1. The command characters in the **pdu** are characters from the standard ASCII character set.

Command character	B	C	D	F	I	R	S	T	W
ASCII	42H	43H	44H	46H	49H	52H	53H	54H	57H

2. The most significant byte is transmitted first and the least significant byte finally when transmitting Flow, Disp Vol, Pause time, Copy No., Back suction time and so on.

3. State1: state byte 1.

Bit 0 – start / stop bit, 1 to start the pump, 0 to stop the pump.

Bit 1 – cw / ccw bit, 1 to run in cw, 0 to run in ccw.

Bit 2 – prime bit, 1 to prime the pump at the max speed 100 rpm.

4. Default **addr**: default pump address (i.e. Pump I.D.#): 1.

5. Pump head No.:

1: YZ2515 2: YZ1515 3: DG (6-roller) 4: DG (10-roller)

6. Pump head No.-Tube No.-Tubeing ID

Pump Head No.	1: YZ1515						
Tube No.	1	2	3	4	5	6	7
Tubing ID (mm)	0.8	1.6	2.4	3.1	4.8	6.4	7.9

Pump Head No.	2: YZ2515						
Tube No.	1	2	3	4			
Tubing ID (mm)	4.8	6.4	7.9	9.6			

Pump Head No.	3: DG (6-roller)				4: DG (10-roller)				
Tube No.	1	2	3	4	5	6	7	8	9

LONGER RS485 Protocol for BT100-1F

Tubing ID (mm)	0.13	0.25	0.51	1.02	1.65	2.00	2.40	2.79	3.17
----------------	------	------	------	------	------	------	------	------	------

7. Examples

a. Write Dispensing Parameter

Control computer command string:

E9 01 0E 57 44 00 00 03 E8 00 00 C8 05 F5 E1 00 00 0A 24

- The above command string from control computer will set dispensing parameter of pump 1 as follows: set Disp Vol to 10.0 mL, set Copy No. to 200, set Disp Flow to 100.0 mL/min, set Pause time to 1.0 s.

- When transmitting a command string, E8H is replaced by E8H 00H.

Pump response:


E9 01 02 57 44 10

b. Read Running Parameter of Flow Control Mode

Control computer command string:

E9 01 02 52 46 17

Pump response:



52 46 – RF

0E E6 B2 80 – show flow (250000000 nL/min=250.0 mL/min)

02 – stop state, run in cw

c. Writing Pump Head and Tubing

Control computer command string:

E9 01 04 57 54 02 02 06

Pump response:

E9 01 02 57 54 00

- The above command string from control computer will set pump head and tubing of pump 1 as follows: pump head to YZ2515, Tubing I.D. to 6.4mm.

8. DB-15 External Control Interface

Pin 2 - RS485 B

Pin 3 - RS485 A