

Accumax Compatibility Chart vs Accumax Pipettes

Non-Filter tips					Compatible Volumes (µl)	
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10µl			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10-B	AT-10-R	AT-10-K	0.1 - 2.5µl 0.5-10 µl 2-20 µl	2.5 µl	0.5 - 10 µl 2-20 µl
Sterile	AT-10-S-B	AT-10-S-R	AT-10-S-K			
Low retention	AT-10-L-B	AT-10-L-R	AT-10-L-K			
Sterile, Low retention	AT-10-S-L-B	AT-10-S-L-R	AT-10-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XL-B	AT-10XL-R	AT-10XL-K	0.1 - 2.5µl 0.5-10 µl 2-20 µl	2.5 µl	0.5 - 10 µl 2-20 µl
Sterile	AT-10XL-S-B	AT-10XL-S-R	AT-10XL-S-K			
Low retention	AT-10XL-L-B	AT-10XL-L-R	AT-10XL-L-K			
Sterile, Low retention	AT-10XL-S-L-B	AT-10XL-S-L-R	AT-10XL-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XN			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XN-B	AT-10XN-R	AT-10XN-K	0.1 - 2.5µl 0.5-10 µl 2-20 µl	2.5 µl	0.5 - 10 µl 2-20 µl
Sterile	AT-10XN-S-B	AT-10XN-S-R	AT-10XN-S-K			
Low retention	AT-10XN-L-B	AT-10XN-L-R	AT-10XN-L-K			
Sterile, Low retention	AT-10XN-S-L-B	AT-10XN-S-L-R	AT-10XN-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	200µl			Variable	Fixed	Multi
	Catalog #					
Standard	AT-200-B	AT-200-R	AT-200-K	5-50 µl 50-200 µl 10-100 µl 20-200 µl	5 µl, 10 µl, 20 µl, 25 µl, 50 µl, 100 µl, 200 µl	5-50 µl 10-100 µl 20-200 µl
Sterile	AT-200-S-B	AT-200-S-R	AT-200-S-K			
Low retention	AT-200-L-B	AT-200-L-R	AT-200-L-K			
Sterile, Low retention	AT-200-S-L-B	AT-200-S-L-R	AT-200-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	200XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-200XL-B	AT-200XL-R	AT-200XL-K	5-50 µl 50-200 µl 10-100 µl 20-200 µl	5 µl, 10 µl, 20 µl, 25 µl, 50 µl, 100 µl, 200 µl	5-50 µl 10-100 µl 20-200 µl 30-300 µl
Sterile	AT-200XL-S-B	AT-200XL-S-R	AT-200XL-S-K			
Low retention	AT-200XL-L-B	AT-200XL-L-R	AT-200XL-L-K			
Sterile, Low retention	AT-200XL-S-L-B	AT-200XL-S-L-R	AT-200XL-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	1000µl			Variable	Fixed	Multi
	Catalog #					
Standard	AT-1000-B	AT-1000-R	AT-1000-K	200-1000 µl 100-1000 µl	250 µl 500 µl 1000 µl	
Sterile	AT-1000-S-B	AT-1000-S-R	AT-1000-S-K			
Low retention	AT-1000-L-B	AT-1000-L-R	AT-1000-L-K			
Sterile, Low retention	AT-1000-S-L-B	AT-1000-S-L-R	AT-1000-S-L-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	1000XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-1000XL-B	AT-1000XL-R	AT-1000XL-K	200-1000 µl 100-1000 µl	250 µl 500 µl 1000 µl	
Sterile	AT-1000XL-S-B	AT-1000XL-S-R	AT-1000XL-S-K			
Low retention	AT-1000XL-L-B	AT-1000XL-L-R	AT-1000XL-L-K			
Sterile, Low retention	AT-1000XL-S-L-B	AT-1000XL-S-L-R	AT-1000XL-S-L-K			
Filter tips					Compatible Volumes	
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10µl			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10-F-B	AT-10-R	AT-10-F-K	0.1 - 2.5µl 0.5-10 µl	2.5 µl	0.5 - 10 µl
Sterile	AT-10-S-F-B	AT-10-S-R	AT-10-S-F-K			
Low retention	AT-10-L-F-B	AT-10-L-R	AT-10-L-F-K			
Sterile, Low retention	AT-10-S-L-F-B	AT-10-S-L-R	AT-10-S-L-F-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XL			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XL-F-B	AT-10XL-R	AT-10XL-F-K	0.1 - 2.5µl 0.5-10 µl 2-20 µl	2.5 µl	0.5 - 10 µl 2-20 µl
Sterile	AT-10XL-S-F-B	AT-10XL-S-R	AT-10XL-S-F-K			
Low retention	AT-10XL-L-F-B	AT-10XL-L-R	AT-10XL-L-F-K			
Sterile, Low retention	AT-10XL-S-L-F-B	AT-10XL-S-L-R	AT-10XL-S-L-F-K			
Packing style:	Bulk	Rack	Reloading stack			
Volume Range:	10XN			Variable	Fixed	Multi
	Catalog #					
Standard	AT-10XN-F-B	AT-10XN-R	AT-10XN-F-K	0.1 - 2.5µl 0.5-10 µl	2.5 µl	0.5 - 10 µl
Sterile	AT-10XN-S-F-B	AT-10XN-S-R	AT-10XN-S-F-K			
Low retention	AT-10XN-L-F-B	AT-10XN-L-R	AT-10XN-L-F-K			
Sterile, Low retention	AT-10XN-S-L-F-B	AT-10XN-S-L-R	AT-10XN-S-L-F-K			

Volume Range:	20µl			Variable
	Catalog #			
Standard	AT-20-F-B	AT-20-F-R	AT-20-F-K	
Sterile	AT-20-S-F-B	AT-20-S-F-R	AT-20-S-F-K	
Low retention	AT-20-L-F-B	AT-20-L-F-R	AT-20-L-F-K	
Sterile, Low retention	AT-20-S-L-F-B	AT-20-S-L-F-R	AT-20-S-L-F-K	

Fixed
5 µl 10 µl 20 µl

Multi

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	50µl			Variable
	Catalog #			
Standard	AT-50-F-B	AT-50-F-R	AT-50-F-K	5-50 µl
Sterile	AT-50-S-F-B	AT-50-S-F-R	AT-50-S-F-K	
Low retention	AT-50-L-F-B	AT-50-L-F-R	AT-50-L-F-K	
Sterile, Low retention	AT-50-S-L-F-B	AT-50-S-L-F-R	AT-50-S-L-F-K	

Fixed
25 µl 50 µl

Multi
5-50 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	100µl			Variable
	Catalog #			
Standard	AT-100-F-B	AT-100-F-R	AT-100-F-K	10-100 µl
Sterile	AT-100-S-F-B	AT-100-S-F-R	AT-100-S-F-K	
Low retention	AT-100-L-F-B	AT-100-L-F-R	AT-100-L-F-K	
Sterile, Low retention	AT-100-S-L-F-B	AT-100-S-L-F-R	AT-100-S-L-F-K	

Fixed
100 µl

Multi
10 - 100 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	200µl			Variable
	Catalog #			
Standard	AT-200-F-B	AT-200-F-R	AT-200-F-K	50-200 µl
Sterile	AT-200-S-F-B	AT-200-S-F-R	AT-200-S-F-K	
Low retention	AT-200-L-F-B	AT-200-L-F-R	AT-200-L-F-K	
Sterile, Low retention	AT-200-S-L-F-B	AT-200-S-L-F-R	AT-200-S-L-F-K	

Fixed
200 µl

Multi
20 - 200 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	200XL			Variable
	Catalog #			
Standard	AT-200XL-F-B	AT-200XL-F-R	AT-200XL-F-K	50-200 µl
Sterile	AT-200XL-S-F-B	AT-200XL-S-F-R	AT-200XL-S-F-K	
Low retention	AT-200XL-L-F-B	AT-200XL-L-F-R	AT-200XL-L-F-K	
Sterile, Low retention	AT-200XL-S-L-F-B	AT-200XL-S-L-F-R	AT-200XL-S-L-F-K	

Fixed
200 µl

Multi
20 - 200 µl 30-300 µl

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	1000µl			Variable
	Catalog #			
Standard	AT-1000-F-B	AT-1000-F-R	AT-1000-F-K	200-1000 µl 100-1000 µl
Sterile	AT-1000-S-F-B	AT-1000-S-F-R	AT-1000-S-F-K	
Low retention	AT-1000-L-F-B	AT-1000-L-F-R	AT-1000-L-F-K	
Sterile, Low retention	AT-1000-S-L-F-B	AT-1000-S-L-F-R	AT-1000-S-L-F-K	

Fixed
250 µl 500 µl 1000 µl

Multi

Packing style:	Bulk	Rack	Reloading stack	
Volume Range:	1000XL			Variable
	Catalog #			
Standard	AT-1000XL-F-B	AT-1000XL-F-R	AT-1000XL-F-K	200-1000 µl 100-1000 µl
Sterile	AT-1000XL-S-F-B	AT-1000XL-S-F-R	AT-1000XL-S-F-K	
Low retention	AT-1000XL-L-F-B	AT-1000XL-L-F-R	AT-1000XL-L-F-K	
Sterile, Low retention	AT-1000XL-S-L-F-B	AT-1000XL-S-L-F-R	AT-1000XL-S-L-F-K	

Fixed
250 µl 500 µl 1000 µl

Multi