

# Innovalene MED 55

## MEDICAL PUMP TUBING

Lower Shore A 55 durometer tubing is specially formulated for medical peristaltic pump use. Ultra smooth inner surface effectively reduces unwanted fluid residue.

Medical



### Features

Transparent  
Ideal for contact with bodily fluids  
Non DEHP  
Peristaltic pump use  
Non-wetting surface permits complete drainage  
Shore A 55 durometer  
Customized size and length

### Regulations

ISO10993  
RoHS

### Applications

IVD  
Blood and IV solutions  
Dialysis equipment  
Inhalation equipment  
Chemotherapy drug delivery  
Peristaltic pumps  
Kidney dialysis pumps



### Innovalene MED 55

Non DEHP tubing has been formulated specifically to meet the demanding challenges of applications such as dialysis, chemotherapy drug delivery, minimally invasive surgeries and peristaltic pump use in analytical and medical devices.



+86 (0)512 58637230  
+86 17317862507



info@innovapure.com  
www.innovapure.com



Innovation Park, Leyu town  
Zhangjiagang, Suzhou, China

## Innovalene™ MED 55 inventoried sizes

Model #	Tube #	I.D. (mm)	I.D. (Inch)	O.D. (mm)	O.D. (Inch)	Wall thickness (mm)	Wall thickness (Inch)	Length / Box* (M)	Length / Box* (Feet)
L305-201-0159-0476	14#	1.59	1/16	4.76	3/16	1.59	1/16	15.2	50
L305-201-0238-0556	19#	2.38	3/32	5.56	7/32	1.59	1/16	15.2	50
L305-201-0318-0635	16#	3.18	1/8	6.35	1/4	1.59	1/16	15.2	50
L305-201-0318-0794	120#	3.18	1/8	7.94	5/16	2.38	3/32	15.2	50
L305-201-0476-0794	25#	4.76	3/16	7.94	5/16	1.59	1/16	15.2	50
L305-201-0476-0953	15#	4.76	3/16	9.53	3/8	2.38	3/32	15.2	50
L305-201-0635-0953	17#	6.35	1/4	9.53	3/8	1.59	1/16	15.2	50
L305-201-0635-1111	24#	6.35	1/4	11.11	7/16	2.38	3/32	15.2	50
L305-201-0794-1111	18#	7.94	5/16	11.11	7/16	1.59	1/16	15.2	50
L305-201-0794-1270	35#	7.94	5/16	12.70	1/2	2.38	3/32	15.2	50
L305-201-0953-1270	96#	9.53	3/8	12.70	1/2	1.59	1/16	15.2	50
L305-201-0953-1429	36#	9.53	3/8	14.29	9/16	2.38	3/32	15.2	50
L305-201-0953-1588	73#	9.53	3/8	15.88	5/8	3.18	1/8	15.2	50
L305-201-1111-1429		11.11	7/16	14.29	9/16	1.59	1/16	15.2	50
L305-201-1270-1905	82#	12.70	1/2	19.05	3/4	3.18	1/8	15.2	50

\* Other sizes and length can be customized

## Innovalene™ MED 55 physical properties

Physical Properties	Test Method	Unit	Typical Value
Color	-		Transluscent
Shore A, 5 sec.	ASTM D2240		55
Specific gravity	ASTM D792	g/cm <sup>3</sup>	1.16
Tensile strength	ASTM D638	Mpa	10.00
Elongation at 100%	ASTM D638	%	400
Tear strength	ASTMD1004-94	kN/m	22
Compression set, 23C 22hrs	ASTM D395-98-B	%	25

Innovalene® MED 55 chemical compatibility chart is available on [www.innovapure.com](http://www.innovapure.com) or call Innovapure sales representative.



### Important

It is the user's responsibility to ensure the suitability and safety of Suzhou Innovapure Co., Ltd. products for all intended uses. Laboratory and clinical tests must be conducted in accordance with applicable regulatory requirements in order to determine the safety and effectiveness for use of tubing in any particular application.